






# Biosearch

## Customer Bottle Catalogue

Picture	Product Code	Description	Use	Price/Quantity
	D1000	<ul style="list-style-type: none"> <li>•1 Litre PET Bottle</li> <li>•Blue Cap</li> <li>•Dosed with 180mg/L of Sodium Thiosulphate</li> <li>•Irradiated</li> </ul>	legionella	£28.80 Pack Size: 48
	D500	<ul style="list-style-type: none"> <li>•500ml PET Bottle</li> <li>•Blue Cap</li> <li>•Dosed with 180mg/L of Sodium Thiosulphate</li> <li>•Irradiated</li> </ul>	bacti+ pseud	£63.00 Pack Size: 126
	D250	<ul style="list-style-type: none"> <li>•250ml PET Bottle</li> <li>•Blue Cap</li> <li>•Dosed with 180mg/L of Sodium Thiosulphate</li> <li>•Irradiated</li> </ul>	bacti	£84.00 Pack Size: 210
	D1000EDTA	<ul style="list-style-type: none"> <li>•1 Litre PET Bottle</li> <li>•Red Cap</li> <li>•Dosed with &gt;10mg/l of EDTA</li> <li>•Irradiated</li> </ul>	legionella in Cu/Hg	£28.80 Pack Size: 48
	D500EDTA	<ul style="list-style-type: none"> <li>•500ml PET Bottle</li> <li>•Red Cap</li> <li>•Dosed with &gt;10mg/l of EDTA</li> <li>•Irradiated</li> </ul>	Bacti+Pseud in Cu/Hg	£63.00 Pack Size: 126
	D250EDTA	<ul style="list-style-type: none"> <li>•250ml PET Bottle</li> <li>•Red Cap</li> <li>•Dosed with &gt;10mg/l of EDTA</li> <li>•Irradiated</li> </ul>	Bacti in Cu/Hg	£84.00 Pack Size: 210

# Biosearch

## Bottle - Catalogue

Picture	Product Code	Description	Use	Price/Quantity
	PET1000	<ul style="list-style-type: none"> <li>•1 Litre PET Bottle</li> <li>•Natural Cap</li> <li>•Un-Dosed</li> </ul>	Chem- All	£19.20 Pack Size: 48
	PET500	<ul style="list-style-type: none"> <li>•500ml PET Bottle</li> <li>•Natural Cap</li> <li>•Un-Dosed</li> </ul>	Chem- pH, Turb, Cond, Anions, TSS, Hardness, TOC, Alk	£28.80 Pack Size: 96
	BIOS 0203SQ	<ul style="list-style-type: none"> <li>•250ml PET Bottle</li> <li>•Natural Cap</li> <li>•Un-Dosed</li> </ul>	Chem- pH, Turb, Cond, Anions	£63.00 Pack Size: 210
	CP125	<ul style="list-style-type: none"> <li>•125ml HDPE Bottle</li> <li>•Natural Cap</li> <li>•Un-Dosed</li> </ul>	Chem- Metals	£105.00 Pack Size: 350
	U500	<ul style="list-style-type: none"> <li>•500ml PET Bottle</li> <li>•Green Cap</li> <li>•Un-Dosed</li> <li>•Irradiated</li> </ul>	Bacti from Ultrapure Process Systems	£50.40 Pack Size: 126